

GTP-Binding Protein Fragment, Gs alpha

Chemical Properties

CAS No. :


Formula: C85H144N30O24S

Molecular Weight: 2002.3

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

GTP-Binding Protein Fragment, Gs alpha is a polypeptide molecule with the sequence CKQLQKDKQVYRATHR.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.4994 mL	2.4971 mL	4.9943 mL
5 mM	0.0999 mL	0.4994 mL	0.9989 mL
10 mM	0.0499 mL	0.2497 mL	0.4994 mL
50 mM	0.010 mL	0.0499 mL	0.0999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481