

## Orexin B, canine

## Chemical Properties

CAS No. :


Formula: C125H214N44O34S

Molecular Weight: 2909.38

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Orexin B, canine is a polypeptide molecule with the sequence  
RPGPPGLQGRLQRLQASGNHAAGILTM-NH<sub>2</sub>.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3437 mL	1.7186 mL	3.4372 mL
5 mM	0.0687 mL	0.3437 mL	0.6874 mL
10 mM	0.0344 mL	0.1719 mL	0.3437 mL
50 mM	0.0069 mL	0.0344 mL	0.0687 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481