

MC-Val-Lys(mPEG23)-PAB-MMAE

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

MC-Val-Lys(mPEG23)-PAB-MMAE is a part of the antibody-conjugated active molecule ADC, featuring degradable and hydrophilic modification characteristics, and can be used for the synthesis of ADC compounds.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481