

DBCO-PEG4-GGFG-DX8951

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

DBCO-PEG4-GGFG-DX8951 is a click-reactive Drug-linker that contains a DBCO group and can undergo a strain-promoted acetylene-nitride cycloaddition reaction (SPAAC) with molecules containing an Azide group. DX8951 is a DNA topoisomerase I (topoisomerase I) inhibitor. DBCO-PEG4-GGFG-DX8951 can be used in the synthesis of antibody-drug conjugates (ADCs).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481