

SPDB-DM1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	SPDB-DM4 is a part of the ADC. It is composed of DM1 (a microtubule inhibitor, also an antibody-encodable aristolochic acid alkaloid, used to overcome the systemic toxicity associated with aristolochic acid and enhance tumor-specific delivery) and SPDB, and exhibits highly efficient anti-tumor activity.
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481