

C-lock-GGFG-Exatecan

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

C-lock-GGFG-Exatecan is a thiol-reactive drug linker. GGFG is a peptidase-activatable peptide. After entering the lesion, it can enable the release and activation of the payload Exatecan. C-lock-GGFG-Exatecan is commonly used in the synthesis of antibody-drug conjugates (ADCs).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481