

MC-VA-PAB-MMAD

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

MC-VA-PAB-MMAD is a thiol-reactive Drug-linker. MMAD is an effective microtubule (tubulin) inhibitor, which is commonly used in the synthesis of antibody-drug conjugates (ADCs).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481