

MPB-VC-PAB-DM1

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

MPB-VC-PAB-DM1 is a thiol-reactive Drug-linker. DM1 is a microtubule inhibitor, which is commonly used in the synthesis of antibody-drug conjugates (ADCs).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481