

MA-PEG4-VC-PAB-DMEA-duocarmycin TM

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

MA-PEG4-VC-PAB-DMEA-duocarmycin TM is a thiol-reactive Drug-linker. Duocarmycin TM is an effective antibiotic with anti-tumor activity, and is commonly used in the synthesis of antibody-drug conjugates (ADC).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481