

## DBCO-PEG4-VC-PAB-DMEA-Duocarmycin SA

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

DBCO-PEG4-VC-PAB-DMEA-Duocarmycin SA is a click-reactive Drug-linker that contains a DBCO group and can undergo strain-promoted acylnitridation addition reaction (SPAAC) with molecules containing azide groups. Duocarmycin SA is a highly effective and orally active cytotoxic agent. DBCO-PEG4-VC-PAB-DMEA-Duocarmycin SA can be used in the synthesis of antibody-drug conjugates (ADC).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481