

## Gly-Gly-Gly-MMAF

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Gly-Gly-Gly-MMAF is a transpeptidase reaction-type Drug-linker, featuring an innovative highly stable and degradable linker. MMAF is an effective tubulin polymerization inhibitor, commonly used in the synthesis of antibody-drug conjugates (ADCs).

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481