

NH2-VC-PAB-DMEA-PNU159682

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

NH2-VC-PAB-DMEA-PNU159682 is a transpeptidase reaction-type Drug-linker with an innovative highly stable and degradable linker. PNU159682 is a DNA topoisomerase II (Topo II) inhibitor with excellent cytotoxicity, which can be used in the synthesis of antibody-drug conjugates (ADCs).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481