

PNU159682-EDA-D4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

PNU159682-EDA-D4 is a deuterium-labeled drug linker used for antibody-drug conjugates (ADCs). It is formed by combining PNU159682-D4 (a potent DNA topoisomerase II inhibitor) with EDA. It is commonly used in the synthesis of ADC compounds. The deuterium labeling can be used as a quantitative tracer during the drug development process.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481