

## CD20 Protein-VLP, Rat, Recombinant (His)

### General Information

Synonyms: Ms4a1

Protein Construction: 1-298 aa

Species: Rat

Expression Host: HEK293 Cells

Accession: D4A4Y3

Molecular Weight: 34.1 kDa (Predicted)

AA Sequence: MTTPRNSVSGPFPTEPTKGPLAMQPAQKVIPRRPSSLVGPTQSFFMRESKALGAVQIMNGLFHISLGLLMIPTGVFAPICLSVWYPLWGGIMYIISGSLLA AAAAEKTSRKSLVKAKVIMSSLSLFAAISGIILSIMDILNITVSHFLKMERLELIKTPKLYVDIYNCEPSNSSEKNSPSTQYCNSVQSVFLGILSVMLISAFFQKLVTAGVVENEWKRVCSRPKSNVLLSAGEKKEQTIKMKEEIIELSGVSSQPKNEEEIEIIPVQEEEEEEAEINFPAPPQE QESLPVENEISP

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: The purity information is not available for VLPs proteins.

Endotoxin: Not tested.

Formulation: Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 55 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481