

## SORD Protein, Mouse, Recombinant (His)

### General Information

Synonyms:	Sdh1;Xylitol dehydrogenase (XDH);Sorbitol dehydrogenase;Polyol dehydrogenase;SDH;SORD;L-iditol 2-dehydrogenase;Sord
Protein Construction:	2-357 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q64442
Molecular Weight:	45.0 kDa (Predicted)
AA Sequence:	AAPAKGENLSLVVHGPDIRLENYPIPELGPNDVLLKMHSVIGICSDVHYWEHGRIGDFVVKKPMVLGHEAA GTVTKVGELVKHLKPGDRVAIEPGVPREVDEYCKIGRYNLTPTIFFCATPPDDGNLCRFYKHNADFCYKLPDSV TFEEGALIEPLSVGIYACRRGSVSLGNKVLVCGAGPVMVTLVAKAMGAAQVVVTDLSASRLTKAKEVGADF TIQVGKETPQEIASKVESLLGSKPEVTIECTGAESSVQTGIYATHSGGTLVIVGMGAEMVNLPLVHAAIREVDIK GVFRYCNTWPMASMLASKTLNVKPLVTHRFPLEKAVEAFETAKKGVGLKVMIKCDPNDQNP

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

Reconstitution:	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Stability & Storage:	Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 107 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>
Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481