

## IA2 Protein, Mouse, Recombinant (His)

### General Information

Synonyms: PTP IA-2; Ptpn; R-PTP-N; Ptp35; Receptor-type tyrosine-protein phosphatase-like N

Protein Construction: 601-979 aa

Species: Mouse

Expression Host: E. coli

Accession: Q60673

Molecular Weight: 49.9 kDa (Predicted)

AA Sequence:

RHHSRQRDKERLAALGPEGAHGDTTFEYQDLCRQHMTKSLFNRAEGQPEPSRVSSVSSQFSDAAQASPSS  
HSSTPSWCEEPAQANMDISTGHMILAYMEDHLRNRDLAKEWQALCAYQAEPNTCAAAQDESNIKKNRHP  
DFLPYDHARIKLVESPSRSDYINASPIIEHDPMPAYIATQGPLSHTIADFWQMVWESGCTVIVMLTPLVED  
GVKQCDRYWPDEGSSLYHVYEVNLVSEHIWCEDFLVRSFYLNKLNLTQETRTLQTFHFLSWPAEGTPASTRPLL  
DFRRKVNKCYRGRSCPIIVHCSDGAGRTGTYLIDMVLNRMAGVKEIDIAATLEHVRDQRPGLVRSKDQFEFA  
LTAVAEVNAALKALPQ

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: Not tested.

Formulation: Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 110 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481