

## AKR1A1 Protein, Mouse, Recombinant (His & Myc)

### General Information

Synonyms:	Glucuronolactone reductase;Aldehyde reductase;Akr1a1;Akr1a4;Alcohol dehydrogenase [NADP(+)];S-nitroso-CoA reductase (ScorR);Glucuronate reductase;Aldo-keto reductase family 1 member A1
Protein Construction:	2-325 aa
Species:	Mouse
Expression Host:	E. coli
Accession:	Q9JII6
Molecular Weight:	43.9 kDa (Predicted)
AA Sequence:	TASSVLLHTGQKMPLIGLGTWKSEPGQVKA AIKHALSAGYRHIDCASVYGNETEIGEALKESVSGKAVPREEL FVTSKLWNTKHHPEDVEPALRRTLADLQLEYLDLYLMHWPYAFERGDNPFKKNADGTVRYDSTHYKETWKA LEVLVAKGLVKALGLSNFNSRQIDDVLSVASVRPAVLQVECHPYLAQNELIAHCHARGLEVTAYSPLGSSDRA WRHPDEPVLLEEPVVLALAEKHGRSPAQILLRWQVQRKVICIPKSINPSRILQNIQVDFDFTFSPEEMKQLDALN KNWRYIVPMITVDGKRVPRDAGHPLYPFNDPY

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 149 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481