

## EPOR Protein, Macaca mulatta, Recombinant (His)

### General Information

Synonyms: Erythropoietin receptor; EPOR; EPO-R

Protein Construction: 26-250 aa

Species: Rhesus macaque

Expression Host: P. pastoris (Yeast)

Accession: F7GNY3

Molecular Weight: 26.9 kDa (Predicted)

AA Sequence:

PPPNLPDPKFESKAALLVARGPEELLCFTEILEDLVCFWEEAASAGVSPDNYSFSYHLEDEPWKLCRLHQAPT  
ARGAVRFWCSLPTADTSSFVPLELRVTAASGAPRYHRIIHINEVLLDAPVGLVARLADEGGHIVLRWLPPPEA  
PMTSHIRYEVDVSAGNGAGSVQRVEILEGRTECVLSNLRGRTRYTF AVRARMAEPSFGGFWSAWSEPVSLTLP  
SDLDP

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 95% as determined by SDS-PAGE.

Endotoxin: Not tested.

Formulation: Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 160 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481