

## PTH1R Protein, Macaca fascicularis, Recombinant (His)

## General Information

Synonyms:	PTH1R;Parathyroid hormone/parathyroid hormone-related peptide receptor;PTH/PTHrP type I receptor;Parathyroid hormone 1 receptor
Protein Construction:	29-188 aa
Species:	Cynomolgus monkey
Expression Host:	HEK293 Cells
Accession:	A0A2K5V724
Molecular Weight:	19.9 kDa (Predicted)
AA Sequence:	DDVMTKEEQIFLLHRAQAQCEKRLKEVLQRPANIMESDKGWTSTSTSGKPRKDKPSGKLYPESEEDKEAPTGS RYRGRPCLPEWDHILCWPLGAPGEVAVPCPDYIYDFNHKGHAYRRCDRNGSWELVPGHNRTWANYSECV KFLTNETREREVFDRLG

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 165 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481