

SLFN11 Protein, Human, Recombinant (aa 1-901, His)

General Information

Synonyms: Schlafen family member 11;SLFN11

Protein Construction: 1-901 aa

Species: Human

Expression Host: E. coli

Accession: Q7Z7L1

Molecular Weight: 106.9 kDa (Predicted)

AA Sequence:

MEANQCPLVVEPSYPDLVINVGEVTLGEENRKKLQKIQRDQEKERMRAACALLNSGGGVIRMAKKVEHPVE
MGLDLEQSLRELIQSSDLQAFFETKQQGRCFYIFVKSWSGGPFPEDRSVKPRCLSSSLYRRSETSVRSMDSR
EAFCLKTKRKP KILEEGPFHKIHKGVYQELPNSDPADPNSDPADLIFQKDYLEYGEILPFPESQLVEFKQFSTK
HFQEYVKRTIPEYVPAFANTGGGYLFIGVDDKSREVLGCAKENVDPDSLRRKIEQAIYKLPVHFCQPQRPIFTF
LKIVNVLKRGEYGYACMIRVNPFCFAVSEAPNSWIVEDKYVCSLTTEKVVGMMDTDPDLLQLSEDFECQL
SLSSGPPLSRPVYSKKGLEHKKELQQLFSVPPGYLRYTPESLWRDLISEHRGLEELINKQMOPFFRGILIFSRS
WAVDLNLQEKPGVICDALLIAQNSTPILYTLIREQDAEGQDYCTRFTAFTLKQKLVNMGGYTGKVCVRKVLCL
SPESSAEALEAAVSPMDYPASYSLAGTQHMEALLQSLVIVLLGFRSLLSDQLGCEVLNLLTAQQYEIFSRSLRK
NRELFVHGLPGSGKTIMAMKIMEKIRNVFHCEAHRILYVCENQPLRNFISDRNICRAETRKTFLRENFEHIQHIVI
DEAQNFRTEGDWYGKAKSITRRAKGGPGILWIFLDYFQTSHLDCSGLPPLSDQYPREELTRIVRNADPIAKYL
QKEMQVIRSNPSFNIPITGCELVFPEAEWSQGVQGLTRIKKYLTVEQIMTCVADTCRRFFDRGYSPKDVAVLVST
AKEVEHYKYELLKAMRKKRVVQLSDACDMLGDHIVLDSVRRFSGLERSIVFGIHPRTADPAILPNVLICLASRA
KQHLYIFPWGGH

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: Not tested.

Formulation: Lyophilized from PBS, 6% Trehalose, pH 7.4

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 199 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481