

## NLRP3 Protein, Human, Recombinant (aa 131-704, His)

## General Information

Synonyms:	Angiotensin/vasopressin receptor AII/AVP-like;PYRIN-containing APAF1-like protein 1; C1orf7;PYPAF1;NACHT, LRR and PYD domains-containing protein 3;CIAS1;Caterpillar protein 1.1 (CLR1.1);Cold-induced autoinflammatory syndrome 1 protein;NALP3;NLRP3;Cryopyrin
Protein Construction:	131-704 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q96P20
Molecular Weight:	73.9 kDa (Predicted)
AA Sequence:	KMKKDYRKKYRKYVRSRFQCIEDRNARLGESVSLNKRYTRLRLIKEHRSQQEREQELLAIGKTKTCESPVSPIK MELLDPPDDEHSEPVHTVVFQGAAGIGKTLARKMMLDWASGTLYQDRFDYLFYIHCREVSLVTQRSLGDLIM SCCPDPNPIHKIVRKRPSRILFLMDGFDELQGAFFEHLGGLCTDWQKAERGDILLSSLRKLLPEASLLITRPV ALEKLQHLLDHPRHVEILGFSEAKRKEYFFKYFSDEAQAARAFSLIQENEVLFTMCFIPLVCWIVCTGLKQQME SGKSLAQTSKTTTAVYVFFLSSLLQPRGGSQEHLCAHLWGLCSLAADGIWNQKILFEESDLRNHGLQKADV SAFLRMNLFQKEVDCEKFYSFIHMTFQEFFAAMYLLLEEEKEGRTNVPGSRLKLPDRDVTVLENYGKFEKGYL IFVVRFLFGLVNQERTSYLEKKLSCKISQQIRLELLKWIEVKAKAKKLIQPSQLELFYCLYEMQEEDFVQRAMD YFPKIEINLSTRMDHMVSSFCIENCHRVESLSLGLHNPKEEEEEKEGRHDMV

## QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

## Preparation and Storage

## Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

## Stability &amp; Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 234 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

## Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481