

HIF1AN Protein, Human, Recombinant (GST)

General Information

Synonyms:	FIH1;Hypoxia-inducible factor asparagine hydroxylase;Factor inhibiting HIF-1 (FIH-1); Hypoxia-inducible factor 1-alpha inhibitor;HIF1AN
Protein Construction:	2-349 aa
Species:	Human
Expression Host:	P. pastoris (Yeast)
Accession:	Q9NWT6
Molecular Weight:	65.8 kDa (Predicted)
AA Sequence:	AATAAEAVASGSGEPREEAGALGPAWDESQLRSYSFPTRPIRLSQSDPRAEELIENEPPVLTDTNLVYPALK WDLEYLQENIGNGDFSVYSASTHKFLYYDEKKMANFQNFKPRSNREEMKFHEFVEKLQDIQQRGGEERLYLQ QTLNDDTVGRKIVMDFLGFNWNWINKQQGKRGWGQLTSNLLLIGMEGNVTPAHYDEQQNFFAQIKGYKRCIL FPPDQFECLYPYPVHHPCDRQSQVDFDNPDYERFPNFQNVVGYETVVGPGDVLVYIPMYWWHHIESLLNGGI TITVNFWYKGAPTPKRIEYPLKAHQKVAIMRNIEMKMLGEALGNPQEVGPLLNTMIKGRYN

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

Preparation and Storage

Reconstitution:	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Stability & Storage:	Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 267 months. Please avoid multiple freeze-thaw cycles and store products in aliquots. <small>Actual storage temperature shall be subject to the COA.</small>

Shipping:	In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.
-----------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481