

## FGF-5 Protein, Human, Recombinant (His)

### General Information

Synonyms:	Smag-82;Heparin-binding growth factor 5 (HBGF-5);Fibroblast growth factor 5;FGF5;FGF-5
Protein Construction:	21-268 aa
Species:	Human
Expression Host:	E. coli
Accession:	P12034
Molecular Weight:	34.2 kDa (Predicted)
AA Sequence:	HGEKRLAPKGQPGPAATDRNPRGSSSRQSSSSAMSSSSASSSPAASLGSGGLEQSSFQWSPSGRRTGSL YCRVGIGFHLQIYPDGKVNGSHEANMLSVLEIFAVSQGIVGIRGVFSNKFLAMSKKGLHASAKFTDDCKFRER FQENSYNTYASAIHRTEKTGREWYVALNKRKAKRGCSPRVKPKQHISTHFLPRFKQSEQPELSFTVTVPEKKK PPSPIKPKIPLSAPRKNTNSVKYRLKFRFG

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 283 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481