

## HUWE1 Protein, Human, Recombinant (His & Myc)

### General Information

Synonyms:	HECT-type E3 ubiquitin transferase HUWE1; Large structure of UREB1 (LASU1); UREB1; HUWE1; Mcl-1 ubiquitin ligase E3 (Mule); Upstream regulatory element-binding protein 1 (URE-B1; URE-binding protein 1); E3 ubiquitin-protein ligase HUWE1; KIAA0312; ARF-binding protein 1 (ARF-BP1); KIAA1578; Homologous to E6AP carboxyl terminus homologous protein 9 (HectH9); HECT, UBA and WWE domain-containing protein 1
Protein Construction:	4005-4374 aa
Species:	Human
Expression Host:	E. coli
Accession:	Q7Z6Z7
Molecular Weight:	50.7 kDa (Predicted)
AA Sequence:	LERLDEGLRKEDMAVHVRRDHVFEDSYRELHRKSPEEMKNRLYIVFEGEEGQDAGGLLREWYMIISREMFNP MYALFRTSPGDRVITYTINPSSHCNPNHLSYFKFVGRIVAKAVYDNRLLCEYFTRSFYKHLGKSVRYTDMESD YHFYQGLVYLLENDVSTLGYDLTFSTEVQEFVCEVRDLKPNGANILVTEENKKEYVHLVCQMRMTGAIRKQL AAFLEGFYEIIPKRLISIFTEQELELLISGLPTIDDLKSNTTEYHKYQSNSIQIWFWRALRSFDQADRKFLQFV TGTSKVPLQGFAALEGMNGIQKFQIHRDDRSTDRLPSAHTCFNQLDLPAYESFEKLRHMLLLAIQECSEGFGL A

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 90% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 296 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481