

APPL1 Protein, Human, Recombinant (GST)

General Information

| | |
|-----------------------|---|
| Synonyms: | Adapter protein containing PH domain, PTB domain and leucine zipper motif 1;KIAA1428; DIP13A;Dip13-alpha;DCC-interacting protein 13-alpha;APPL1;APPL |
| Protein Construction: | 5-377 aa |
| Species: | Human |
| Expression Host: | E. coli |
| Accession: | Q9UKG1 |
| Molecular Weight: | 70.2 kDa (Predicted) |
| AA Sequence: | DKLPIEETLEDSPQTRSLLGVFEEDATAISNYMNQLYQAMHRIYDAQNELSAATHLTSKLLKEYEKQRFPLGG DDEVMSSTLQQFSKVIDELSSCHAVLSTQLADAMMFPITQFKERDLKEITLKEVFQIASNDHDAAINRYSRLS KKRENDKVKYEVTEDEVYTSRKKQHQTMMHYFCALNTLQYKKKIALLEPLLGYMQAQISFFKMGSENLNEQLE EFLANIGTSVQNVRREMSDIETMQQTIEDLEVASDPLYVDPDPDKFPVNRNLTRKAGYLNARNKTGLVSST WDRQFYFTQGGNLMSQARGDVAGGLAMDIDNCSVMAVDCEDRRYCFQITSFDGKSSILQAESKKDHEEWI CTINNISKQI |

QC Testing

| | |
|----------------------|---|
| Biological Activity: | Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first. |
| Purity: | > 85% as determined by SDS-PAGE. |
| Endotoxin: | Not tested. |
| Formulation: | Lyophilized from PBS, 6% Trehalose, pH 7.4 |

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 307 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481