

## CCNYL1 Protein, Human, Recombinant (His)

### General Information

Synonyms: Cyclin-Y-like protein 1;CCNYL1

Protein Construction: 1-359 aa

Species: Human

Expression Host: E. coli

Accession: Q8N7R7

Molecular Weight: 41.8 kDa (Predicted)

AA Sequence:

MGNTLTCCVSPNASPKLGRRAGSAELYCASDIYEAVSGDAVAVAPAVVEPAELDFGEGEGHHLQHISDREMP  
EDLALESNPSDHPRASTIFLSKSQTDVREKRKSNHLNHVSPGQLTKKYSSCSTIFLDDSTVSQPNLRRTTVKCVT  
LAIYYHIKNRDANRSLDIFDERSHPLTREKVP EEFYFKHDPEHKFIYRFVRTLFSAAQLTAECAIVTLVYLERLLTYA  
EIDICPTNWKRIVLGAILLASKVWDDQAVWNVDYQCILKDITVEDMNEMERHFLELLQFNINVPASVYAKYYF  
DLRSLADDNNLNFLFAPLSKERAQNLEAISRLCEDKDL CRAAMRRSFSADNFIGIQRSKAILS

### QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 85% as determined by SDS-PAGE.

Endotoxin: Not tested.

Formulation: Lyophilized from PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 312 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

*Actual storage temperature shall be subject to the COA.*

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481