

Aspartate Aminotransferase/GOT2 Protein, Human, Recombinant (His & Myc)

General Information

Synonyms:	Plasma membrane-associated fatty acid-binding protein (FABPpm);KYAT4;Kynurenine aminotransferase IV;Glutamate oxaloacetate transaminase 2;Fatty acid-binding protein (FABP-1);Transaminase A;Kynurenine aminotransferase 4;Kynurenine--oxoglutarate transaminase IV;Aspartate aminotransferase, mitochondrial;mAspAT;Kynurenine--oxoglutarate transaminase 4;GOT2
Protein Construction:	30-430 aa
Species:	Human
Expression Host:	E. coli
Accession:	P00505
Molecular Weight:	52.2 kDa (Predicted)
AA Sequence:	SSWVTHVEMGPPDPILGVTEAFKRDTNSKMMNLGVGAYRDDNGKPYVLPVSRKAEAQIAAKNLDKEYLPIG GLAEFCKASAEALALGENSEVLKSGRFVTVQTISGTGALRIGASFLQRFFKFSRDVFLPKPTWGNHTPIFRDAGM QLQGYRYDPKTCGFDFTGAVEDISKIPEQSVLLHACAHNPTGVDPRPEQWKEIATVVKRNLFAFFDMAY QGASGDGDKDAWAVRHFIEQGINVCLCQSYAKNMGLYGERVGAFTMVCKDADEAKRVESQLKILIRPMYS NPPLNGARIAAAILNTPDLRKQWLQEVKVMADRIIGMRTQLVSNLKEGSTHNWQHITDQIGMFCFTGLKPE QVERLIKEFSIYMTKDGRISVAGVTSSNVGYLAHAIHQVTK

QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 85% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 339 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481