

## IL-18 Protein, Canine, Recombinant (C-His)

### General Information

Synonyms:	IL-18; Interleukin-18; Interferon gamma-inducing factor (IFN-gamma-inducing factor); Interleukin-1 gamma (IL-1 gamma); IL18; IGIF
Protein Construction:	37-193 aa
Species:	Canine
Expression Host:	E. coli
Accession:	Q9XSR0
Molecular Weight:	24.9 kDa (Predicted)
AA Sequence:	YFGKLEPKLSIIRNLNDQVLVNEGNQPVFEDMPDSDCTDNAPHTIFIYMYKDSLTRGLAVTISVKYKTMSTLS CKNKTISFQKMSPDSDINDEGNDIIFQRSVPGHDDKIQFESSLYKGFHLACKKENDLFLKILKDKDENGDKSI MFTVQNKS

### QC Testing

Biological Activity:	Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.
Purity:	> 95% as determined by SDS-PAGE.
Endotoxin:	Not tested.
Formulation:	Lyophilized from PBS, 6% Trehalose, pH 7.4

### Preparation and Storage

#### Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

#### Stability & Storage:

Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 381 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

#### Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

**This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use**

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481