

cyanine7 maleimide

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cyanine7 amine is effective for labeling biomolecules such as proteins, nucleic acids, and cells, facilitating localization, tracking, and detection in various bioanalytical techniques. It provides fluorescent signals for visualization. The excitation/emission wavelengths (Ex/Em) are 750/773 nm.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481