


Fenaminstrobin

Chemical Properties

CAS No. :	366815-39-6
Formula:	C ₂₁ H ₂₁ Cl ₂ N ₃ O ₃
Molecular Weight:	434.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fenaminstrobin is a biochemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3024 mL	11.5122 mL	23.0245 mL
5 mM	0.4605 mL	2.3024 mL	4.6049 mL
10 mM	0.2302 mL	1.1512 mL	2.3024 mL
50 mM	0.046 mL	0.2302 mL	0.4605 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481