


3,5-Dibromopyridine

Chemical Properties

CAS No. :	625-92-3
Formula:	C ₅ H ₃ Br ₂ N
Molecular Weight:	236.892
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3,5-Dibromopyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2214 mL	21.1068 mL	42.2137 mL
5 mM	0.8443 mL	4.2214 mL	8.4427 mL
10 mM	0.4221 mL	2.1107 mL	4.2214 mL
50 mM	0.0844 mL	0.4221 mL	0.8443 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481