


1-Boc-3-piperidone

Chemical Properties

CAS No. :	98977-36-7
Formula:	C ₁₀ H ₁₇ NO ₃
Molecular Weight:	199.2469
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1-Boc-3-piperidone is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0188 mL	25.0941 mL	50.1882 mL
5 mM	1.0038 mL	5.0188 mL	10.0376 mL
10 mM	0.5019 mL	2.5094 mL	5.0188 mL
50 mM	0.1004 mL	0.5019 mL	1.0038 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481