


1-Boc-1,2,3,6-tetrahydropyridine-4-boronic acid pinacol ester

Chemical Properties

CAS No. :	286961-14-6
Formula:	C ₁₆ H ₂₈ BNO ₄
Molecular Weight:	309.209
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	1-Boc-1,2,3,6-tetrahydropyridine-4-boronic acid pinacol ester is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.234 mL	16.1702 mL	32.3405 mL
5 mM	0.6468 mL	3.234 mL	6.4681 mL
10 mM	0.3234 mL	1.617 mL	3.234 mL
50 mM	0.0647 mL	0.3234 mL	0.6468 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481