


2-Fluorophenol

Chemical Properties

CAS No. :	367-12-4
Formula:	C ₆ H ₅ FO
Molecular Weight:	112.10
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Fluorophenol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.9206 mL	44.603 mL	89.2061 mL
5 mM	1.7841 mL	8.9206 mL	17.8412 mL
10 mM	0.8921 mL	4.4603 mL	8.9206 mL
50 mM	0.1784 mL	0.8921 mL	1.7841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481