


4-Chloro-pyridin-2-ylamine

Chemical Properties

CAS No. :	19798-80-2
Formula:	C ₅ H ₅ ClN ₂
Molecular Weight:	128.56
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-Chloro-pyridin-2-ylamine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7785 mL	38.8923 mL	77.7847 mL
5 mM	1.5557 mL	7.7785 mL	15.5569 mL
10 mM	0.7778 mL	3.8892 mL	7.7785 mL
50 mM	0.1556 mL	0.7778 mL	1.5557 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481