


4-(Diphenylamino)benzeneboronic acid

Chemical Properties

CAS No. :	201802-67-7
Formula:	C ₁₈ H ₁₆ BN ₂ O ₂
Molecular Weight:	289.136
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-(Diphenylamino)benzeneboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4585 mL	17.2927 mL	34.5853 mL
5 mM	0.6917 mL	3.4585 mL	6.9171 mL
10 mM	0.3459 mL	1.7293 mL	3.4585 mL
50 mM	0.0692 mL	0.3459 mL	0.6917 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481