


## (2-Bromophenyl)boronic acid

## Chemical Properties

CAS No. :	244205-40-1
Formula:	C <sub>6</sub> H <sub>6</sub> BBrO <sub>2</sub>
Molecular Weight:	200.826
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	(2-Bromophenyl)boronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9793 mL	24.8967 mL	49.7934 mL
5 mM	0.9959 mL	4.9793 mL	9.9587 mL
10 mM	0.4979 mL	2.4897 mL	4.9793 mL
50 mM	0.0996 mL	0.4979 mL	0.9959 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481