


## 3-(N-Boc-amino)phenylboronic acid

## Chemical Properties

CAS No. :	380430-68-2
Formula:	C <sub>11</sub> H <sub>16</sub> BN <sub>1</sub> O <sub>4</sub>
Molecular Weight:	237.06
Storage:	Store at RT Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-(N-Boc-amino)phenylboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2183 mL	21.0917 mL	42.1834 mL
5 mM	0.8437 mL	4.2183 mL	8.4367 mL
10 mM	0.4218 mL	2.1092 mL	4.2183 mL
50 mM	0.0844 mL	0.4218 mL	0.8437 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481