


(3-Fluoro-2-formylphenyl)boronic acid

Chemical Properties

CAS No. :	871126-15-7
Formula:	C ₇ H ₆ BF ₃ O ₃
Molecular Weight:	167.93
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	(3-Fluoro-2-formylphenyl)boronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9549 mL	29.7743 mL	59.5486 mL
5 mM	1.191 mL	5.9549 mL	11.9097 mL
10 mM	0.5955 mL	2.9774 mL	5.9549 mL
50 mM	0.1191 mL	0.5955 mL	1.191 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481