


## N-Boc-4-bromopiperidine

## Chemical Properties

CAS No. :	180695-79-8
Formula:	C <sub>10</sub> H <sub>18</sub> BrNO <sub>2</sub>
Molecular Weight:	264.159
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Boc-4-bromopiperidine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7856 mL	18.9279 mL	37.8558 mL
5 mM	0.7571 mL	3.7856 mL	7.5712 mL
10 mM	0.3786 mL	1.8928 mL	3.7856 mL
50 mM	0.0757 mL	0.3786 mL	0.7571 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481