


2-Bromo-4-(trifluoromethyl)pyridine

Chemical Properties

CAS No. :	175205-81-9
Formula:	C ₆ H ₃ BrF ₃ N
Molecular Weight:	225.994
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Bromo-4-(trifluoromethyl)pyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.425 mL	22.1249 mL	44.2497 mL
5 mM	0.885 mL	4.425 mL	8.8499 mL
10 mM	0.4425 mL	2.2125 mL	4.425 mL
50 mM	0.0885 mL	0.4425 mL	0.885 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481