


(3-Hydroxyphenyl)boronic acid

Chemical Properties

CAS No. :	87199-18-6
Formula:	C ₆ H ₇ BO ₃
Molecular Weight:	137.929
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(3-Hydroxyphenyl)boronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.2501 mL	36.2503 mL	72.5005 mL
5 mM	1.450 mL	7.2501 mL	14.5001 mL
10 mM	0.725 mL	3.625 mL	7.2501 mL
50 mM	0.145 mL	0.725 mL	1.450 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481