


2,6-Dichloro-3-cyanopyridine

Chemical Properties

CAS No. :	40381-90-6
Formula:	C ₆ H ₂ Cl ₂ N ₂
Molecular Weight:	172.999
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,6-Dichloro-3-cyanopyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7803 mL	28.9017 mL	57.8035 mL
5 mM	1.1561 mL	5.7803 mL	11.5607 mL
10 mM	0.578 mL	2.8902 mL	5.7803 mL
50 mM	0.1156 mL	0.578 mL	1.1561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481