


## 3-Carboxyphenylboronic acid pinacol ester

## Chemical Properties

CAS No. :	269409-73-6
Formula:	C <sub>13</sub> H <sub>17</sub> BO <sub>4</sub>
Molecular Weight:	248.083
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Carboxyphenylboronic acid pinacol ester is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.031 mL	20.1548 mL	40.3096 mL
5 mM	0.8062 mL	4.031 mL	8.0619 mL
10 mM	0.4031 mL	2.0155 mL	4.031 mL
50 mM	0.0806 mL	0.4031 mL	0.8062 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481