


## 3-Tolylboronic acid

## Chemical Properties

CAS No. : 17933-03-8  
Formula: C<sub>7</sub>H<sub>9</sub>BO<sub>2</sub>  
Molecular Weight: 135.956  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Tolylboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.3551 mL	36.7755 mL	73.551 mL
5 mM	1.471 mL	7.3551 mL	14.7102 mL
10 mM	0.7355 mL	3.6776 mL	7.3551 mL
50 mM	0.1471 mL	0.7355 mL	1.471 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481