


## 6-Bromoquinoline

## Chemical Properties

CAS No. :	5332-25-2
Formula:	C <sub>9</sub> H <sub>6</sub> BrN
Molecular Weight:	208.055
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	6-Bromoquinoline is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8063 mL	24.0315 mL	48.0631 mL
5 mM	0.9613 mL	4.8063 mL	9.6126 mL
10 mM	0.4806 mL	2.4032 mL	4.8063 mL
50 mM	0.0961 mL	0.4806 mL	0.9613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481