


## 4-Bromo-2-methylpyridine

## Chemical Properties

CAS No. :	22282-99-1
Formula:	C <sub>6</sub> H <sub>6</sub> BrN
Molecular Weight:	172.023
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-Bromo-2-methylpyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8133 mL	29.0664 mL	58.1328 mL
5 mM	1.1627 mL	5.8133 mL	11.6266 mL
10 mM	0.5813 mL	2.9066 mL	5.8133 mL
50 mM	0.1163 mL	0.5813 mL	1.1627 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481