


2-Bromo-3-methoxypyridine

Chemical Properties

CAS No. :	24100-18-3
Formula:	C ₆ H ₆ BrNO
Molecular Weight:	188.022
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Bromo-3-methoxypyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3186 mL	26.5929 mL	53.1858 mL
5 mM	1.0637 mL	5.3186 mL	10.6372 mL
10 mM	0.5319 mL	2.6593 mL	5.3186 mL
50 mM	0.1064 mL	0.5319 mL	1.0637 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481