


## Isoquinolin-4-amine

## Chemical Properties

CAS No. :	23687-25-4
Formula:	C <sub>9</sub> H <sub>8</sub> N <sub>2</sub>
Molecular Weight:	144.1732
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Isoquinolin-4-amine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.9363 mL	34.6813 mL	69.3626 mL
5 mM	1.3873 mL	6.9363 mL	13.8725 mL
10 mM	0.6936 mL	3.4681 mL	6.9363 mL
50 mM	0.1387 mL	0.6936 mL	1.3873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481