


## Fmoc-Hyp-OH

## Chemical Properties

CAS No. :	88050-17-3
Formula:	C <sub>20</sub> H <sub>19</sub> NO <sub>5</sub>
Molecular Weight:	353.3686
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Fmoc-Hyp-OH is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8299 mL	14.1495 mL	28.299 mL
5 mM	0.566 mL	2.8299 mL	5.6598 mL
10 mM	0.283 mL	1.4149 mL	2.8299 mL
50 mM	0.0566 mL	0.283 mL	0.566 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481